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TRAUS

HIGHSPEED (ATN400, CAB10LN, CAB10LK)

Please read the instruction manual carefully before operation and keep it for future reference.

Model	ATN400	CAB10LN	CAB10LK
Picture			
Weight	53g	58g	65g
RPM	400,000 rpm \pm 10% (At the recommended drive air pressure)		
Bur Type	Type 3. by ISO 1797, FG Bur (\varnothing 1.6mm)		
Max. Bur Length	21mm		
Max. Working Part Diameter	\varnothing 2mm		
Recommended Drive Air Pressure	0.30Mpa (3.0kgf/cm ²)		
Operating Drive Air Pressure Range	0.25 ~ 0.35Mpa (2.5 ~ 3.5kgf/cm ²)		
Max.Air Consumption	<66NL/min (0.3Mpa)		

Water Pressure	0.20Mpa (2.0kgf/cm ²)		
Recommended water content of spray	>50ml/min		
Chip Air Pressure	0.10 ~ 0.20Mpa (1.0 ~ 2.0kgf/cm ²)		
Optic	-	Glass Rod	Glass Rod
Coupling	CP10NQ (NSK QD-JM4 Type)	CP10LN (NSK PTL-CL-LED Type)	CP10LK (KAVO MULTIflex LED Type)

※ The rotation speed varies slightly depending on the direction of rotation and the shape of the hose used.
※ Coupling: Designed to NSK Coupling Type. (ATN400, CAB10LN.)
※ Coupling: Designed to KAVO Coupling Type. (CAB10LK)

1. Intended Use

: This air-powered dental handpiece is intended for removal of carious material, cavities and crown preparations, removal of fillings, processing of tooth, restoration of surfaces and as a surgical tool for third molar removal procedures. It is designed for use by a trained professional in the field of general dentistry.


2. Components

: Highspeed, Repair Kit (Wrench, Cleaning Pin, Brush), Spray Nozzle, Instruction Manual, Spare O-Ring (CAB10LN)

3. How to connect each part


3.1 Highspeed

- Mounting ATN400 / CAB10LN
 - Insert the coupling into the hose and tighten the hose nut.
 - Insert the highspeed into the coupling
 - Connect the highspeed and the coupling by pulling connection ring
 - Grasp the coupling and connect the highspeed and the coupling to the end
 - Check if the highspeed is pulled out or not
- Isolation
 - Pull the connection ring of the coupling and remove the highspeed.(ATN400 / CAB10LN)
 - Remove the highspeed while holding the coupling and pulling the highspeed. (CAB10LK)

3.2 Bur

- Mounting
 - Insert the bur
 - Press the button on the highspeed head and connect the bur to the end.
 - Release the button
 - Make sure that the bur is not pulled out by pulling out
- Isolation
 - Pull out the bur pressing the button, the highspeed head.
 - Be sure to check if the product is cleaned, lubricated and sterilized before using it.
 - Remove the bur and disassemble the highspeed after checking if the operation stopped for sure.
 - Do not use the bur that exceed the maximum length. The maximum length of the bur is 21mm.
 - When using the bur, follow the instructions provided by the bur manufacturer.
 - Before use, check whether it is compatible with the bur and use it


4. Maintenance

4.1 Cleaning

- Rotate the head in a tap water about 10 seconds with maximum speed to remove the remains (blood, physiologic saline Solution, etc.) and rinse.
- Remove contaminants on the highspeed using a nylon brush after soaking for 20 minutes in Enzymatic Detergent.
- Rinse the highspeed using the tap water for 3 minutes to remove Enzymatic Detergent.
- Wash the highspeed for 10 minutes in the ultrasonic cleaner containing Enzymatic Detergent.
- Rinse the highspeed for 3 minutes in the ultrasonic cleaner containing purified water repeatedly subjected 3 times. (Change the purified water each time.)
- Wipe using dust-absorbent after the product soaked in alcohol(70%) for 30 seconds.


- ※ If the spray is not working well, clean the spray port using the cleaning pin as shown below.
- ※ If there is foreign material around the clean head hole and it cannot be removed, use a brush to remove it.


- ! - Do not clean using thinner, benzene, or other solvents.
- Do not forcibly insert the cleaning pin into the spray port or push it tilted. The injection nozzle may be deformed, which may cause malfunction.
- After washing, perform a visual inspection in a place with sufficient light. If contamination is found on the surface, please re-clean according to our cleaning method.

4.2 Lubrication

- Insert the spray nozzle to hose connection part(the end part of highspeed) and spray lubricant 2~3 seconds.

- ATN400: NSK QD-J connection spray Nozzle
- CAB10LN : NSK MachLite connection spray Nozzle
- CAB10LK : KaVo connection spray Nozzle
- If the chuck is contaminated, it may cause problems in tightening. Please lubricate the chuck as shown below once a week.

- ! - Lubricate the product after each use.(Pana Plus K131014)
- When lubrication, please hold the product can go out by the spray pressure.
- The nozzle tip is not included in a set.

4.3 Sterilization


- Dry the product completely after cleaning.
- Dry naturally at room temperature.
- Put the product in a sterile bag and seal tightly.
- Put the sterile bag in the sterilizer and sterilization for 4 minutes with 132°C, and let it dry for 30 minutes.

Sterilization	Temperature	Time	Dry Time	132°C
Moist heat sterilization (pre-vacuum)	132°C	4 minutes	30 minutes	SSS

- ! - Maximum sterilization number is 250 times.
- Make sure to use the sterilization pouch and seal tightly.
- Make sure to dry the applicable products before sterilization.
- Sterilizing the product with blood or other remains may cause the breakage of the product. Therefore, it is necessary to clean the product before sterilization. Supplied in non-sterile condition. Be sure to sterilize before use.

4.4 Replace the air turbine cartridge

- Mount the bur on the chuck.
- Loosen the head cap with the head cap wrench
- Push the bur from the head to remove the cartridge.
- Clean the inside of the head.
- Put the Cartridge slowly into the head. Be careful not to damage the ball bearing.
- Lock the head cap with the head cap wrench.


- ! - Replace the cartridge when you feel abnormal vibration, noise or heat while test or operation.

4.5 O-ring exchange <CAB10LN>


Replace the O-ring if you notice any water coming out from the exhaust. If left un-replaced, water may leak from the coupling connection.

- Be sure to replace as a set with the spare o-ring of the component.
- Remove the handpiece from the coupling.
 - Remove the 'Handle B'.
 - Take out the old O-ring with your fingers and insert the spare O-ring into the groove. Be careful not to damage the O-ring.
 - Tighten the 'Handle B' firmly to the handpiece.
- 1 set of O-ring consists of 5 O-rings.

4.6 Cleaning Glass Rod

If dirt or foreign Material Adheres to the glass rod, wipe it with a cotton swab dipped in sterilized alcohol.

- ! - When cleaning the glass rod, sharp edges may cause scratches on its Surface, Which may reduce the light transmittance.

5. Warranty

SAESHIN guarantees this product for 1 year from invoice date. Product warranty means that SAESHIN has responsibility for defective material or operation failure. The product warranty does not cover the user misuse, wrong installation, inappropriate maintenance and repair, and normal wear of consumable parts such as cartridges and bearings. In order for SAESHIN to confirm if the warranty is valid it will require checking the operating condition, environment information, serial number and stamped invoice. The warranty will be performed as repair or exchange according to manufacturer's analysis of product and judgement.

Warranty Exclusions

Misuse and Improper handling by users.
Dropping the products during installation, movement or use.
If the unit is damaged due to unspecified consumable parts or accessories.
If the product is not repaired by Saeshin, or the repair shop designated by Saeshin.
Normal wear of consumable parts.
Damage and failure due to natural disaster.
If the product is disassembled or modified at will.

[Precautions]

- Use for any other purpose is prohibited.
- The company did not sterilize the product. Please make sure to sterilize before operation. (Please refer to the 4.3 Sterilization)
- Do not drop the product or force with strong external shocks.
- Do not disassemble or modify the product.
- Check for vibration, heat, noise, and other abnormal symptoms before use.
- If the product operates abnormally, stop immediately and contact the seller for a check.
- Do not remove the air hose while the product is rotating.
- Moisture and impurities may cause malfunction of the product. Please clean and dry with air.
- Make patient safety the top priority when use.
- To Keep maintain it clean, always wipe the product with rubbing alcohol before use.
- Improper use may cause damage to the product. Be sure to use it by a skilled technician or specialist.
- The manufacturer is not responsible for accidents caused by use, not in accordance with the instruction manual.
- The user is responsible for the proper use and maintenance of the product.
- This product should not be disposed with general waste.
- Federal law restricts this device to sale by or on the order of a Dentist.

Use conditions	Shipping/Storage conditions
Temperature: 10~35°C	Temperature: -20~60°C
Humidity: 10~80%	Humidity: 10~90%
Atmospheric pressure: 700~1,060hPa	Atmospheric pressure: 500~1,060hPa

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