

traus

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

TRAUS

## Implant Angle Handpiece

Before the operation, read the instruction manual carefully and keep it for future reference.

Name	TRAUS CRB26LX	TRAUS CRB26XX	TRAUS CRB27LX	TRAUS CRB27XX	TRAUS CRB36LN	TRAUS CRB36NN	TRAUS CRB46LN	TRAUS CRB46NN	TRAUS CRB46LNTi	TRAUS CRB46NNTi
Max. speed (rpm)	2,000	2,000	1,250	1,250	2,000	2,000	2,000	2,000	2,000	2,000
Gear ratio	20:1		32:1		20:1					
Weight	82g	84g	82g	84g	82g	84g	85g	85g	85g	85g
Optic / Non-optic	Optic	Non-optic	Optic	Non-optic	optic	Non-optic	Optic	Non-optic	Optic	Non-optic
Bur replacement type	Button type									
Articles	Low speed angle									
Chuck type	Type 1. By ISO 1797									
Bur dimension	CA bur (Ø 2.35)									
Length of shank	9~12 mm									
Max overall length of rotary instrument	30 mm									
Coupling dimension	Standard dimension of dental handpiece (ISO 3964)									

- The maximum speed is different with a output of the control box.  
(Maximum speed of table is based on 40,000.)

### 1. Intended purpose

The handpiece is used for dental implant surgery equipped with bur.

#### A. Implant placement, including

1. Preparation of the osteotomy site
2. Bone contouring, osteoplasty

#### B. Periodontal surgeries

1. Bone contouring & alveoplasty around living teeth
2. Removal of exostosis

#### C. Bone grafting

1. Preparation of the donor site (for e.g. symphysis and ascending rames etc.)
2. Harvesting autogen living bone
3. Sinus elevation & grafting of alveolar sockets

#### D. Removal and sectioning of teeth and teeth bone for e.g. impacted third molars and complicated extractions

### 2. Components

Angle, Spray nozzle, Instruction manual

### 3. Bur replacement

- Button type : After press the button and remove the bur, replace the bur.

※ Please use ISO standard bur and avoid to use of rusty or curved bur because it can cause excessive wear of the bearing.

※ Please comply with instructions of bur.  
(Safety speed and precautions)

### 4. Maintenance of the product


#### 1) Cleaning

- ① Rotate the head in a tap water about 10 seconds with maximum speed to remove the remains (blood, physiologic saline solution, etc.) and rinse.
- ② Disconnect the angle handpiece to the motor.
- ③ Remove contaminants on the angle handpiece using a nylon brush after soaking for 20 minutes in Enzymatic Detergent.
- ④ Rinse the angle handpiece using the tap water for 3 minutes to remove Enzymatic Detergent.
- ⑤ Wash the angle handpiece for 10 minutes in the ultrasonic cleaner containing Enzymatic Detergent.
- ⑥ Rinse the angle handpiece for 3 minutes in the ultrasonic cleaner containing purified water repeatedly subjected 3 times.  
(Change the purified water each time.)
- ⑦ Wipe using dust-absorbent after the product soaked in alcohol(70%) for 30 seconds.

※ The following products can be cleaned in disassembled state.  
- TRAUS CRB46LN, TRAUS CRB46NN, TRAUS CRB46LNTi, TRAUS CRB46NNTi

#### ① Disassemble the product.

- A. Grip the handle tightly and turn the head to the direction of the arrow.
- B. Pull the head out of the handle.
- C. Pull the gear assy out of the handle.


② - Disassemble the product after separating the motor from the handpiece.

- ② Rinse the disassembled parts with a cleansing solution sufficiently.  
- Please do not rub the disassembled parts with a stuff like a brush.  
A breakdown can be caused by a bristle of a brush in the product.

#### ③ Rinse the cleansing solution sufficiently.


#### ④ Reassemble the product.

- A. Insert the gear assy in the right direction as far as it will go.
- B. Please turn the gear to check if the gear assy and the lower joint are engaging.

C. Point the head's mark at the handle's open mark, and then insert as far as it will go.

- ② - Please check if the gear assy and the head gear are engaging.  
A forced insertion can cause a breakdown of the gear.
- D. Turn the head to the lock mark in the direction like the picture.



#### ⑤ Dry all surfaces, and then remove inside moisture by lubricating.

② - Please be sure to lock the head of the product tightly.  
If not, it can cause an accident.

#### 2) Lubrication

: Insert the nozzle to motor assembly part(the end part of angle) and spray lubricant 2~3 seconds.

- If the continuous operation, please put the oil every an hour.
- When lubrication, please hold the handpiece firmly because the handpiece can go out by the spray pressure.

#### 3) Sterilization

- ① Put the angle in a sterile bag and seal tightly.
- ② Put the sterile bag in the sterilizer and sterilization for 4 minutes with 132°C, and let it dry for 30 minutes.

Sterilization	Temperature	Time	Dry
Moist heat sterilization	132°C	4 minutes	30 minutes

※ This product is designed to withstand 250 times of re-sterilization process.

※ This product recommends the moist heat sterilization (autoclave), other sterile methods are not guaranteed the effectiveness.

※ The product is intended to be wrapped in an FDA cleared sterilization pouch/bag prior to sterilization.

#### [Caution for use]

- This product can be used for intended propose only.
- This product is not in sterilization once release by manufacturer.  
Please certainly do user sterilization before using.  
(Refer tto 'Maintenance of product – sterilization')
- If sterilize the product with rest of blood, the product can be damaged.  
So, please fully clean up and lubricate the product.
- Do not use this product for wrong purpose and be aware of the safety.
- Keep the recommended loading & resting time by manufacturer of control box. If not, patient might get burn of damage by overheating.
- After stopping the motor perfectly, disassemble the product from the motor.
- Disassemble or repair the product only the case which is mentioned in instruction manual.
- Do not drop the product. It can be damaged.
- Please insert the product to the end of E-type motor shaft and use it.
- Verify the product's operation, especially vibration, noise, and temperature rise are normal.
- Please contact the dealer/store or manufacturer when problem occurs.
- The manufacturer does not have any responsibility for any accidents that user did not follow the instruction manual.
- User has responsibility for the use and maintenance of the product.
- US federal law restricts this device to sale by or on the order of dentist.
- If you want to get an instruction manual in other languages, please contact the manufacturer.

Operating condition	Transport and storage condition
Temperature : 10~35°C	Temperature : -20~60°C
Humidity : 10~80%	Humidity : 10~90%
Atmospheric pressure : 700~1,060hPa	Atmospheric pressure : 500~1,060hPa

This product is the medical device.

<b>SAESHIN</b> SAESHIN PRECISION CO., LTD. 52, Secheon-ro 1-gil, Dasa-eup, Dalseong-gun, Daegu, Korea Tel. 82-53-987-2341 / Fax 82-53-580-0999 www.saeshin.com / traus@saeshin.com	<b>MD</b> <b>REP</b> MERIDIUS MEDICAL EUROPE LIMITED Unit 3D, North Point House North point Business Park New Mallow Road CORK, T23AT27, Ireland +353 212066448 SOM-S073P-KE  2023-02-01
---	---